



The
Heart
of Infection
Control

AUSTRALIA & NEW ZEALAND

AGED CARE INFECTION CONTROL

PRODUCT CATALOGUE



Care
for Staff &
Residents



Confidence
in Products &
Protocols



Compliance
with Standards

QUICK LINKS

Patient Care & Infection Control In Aged Care page 3

PATIENT CARE page 4

Patient Hand Hygiene page 4

Patient Personal Hygiene page 4

Patient Contenance Management page 5

Patient Personal Protective Equipment (PPE) page 5

Patient Apparel page 5

INFECTION CONTROL page 6

Hand Hygiene page 6

Masks page 8

Protective Eyewear page 12

Test Methods for Face Masks page 13

Gloves page 14

Gowns page 20

Footwear page 24

Ankle Tops page 24

Knee Highs page 24

Protection for Surfaces & Equipment page 25

Surface Disinfection page 25

Surface Protection page 26

Patient Linen page 29

PATIENT CARE & INFECTION CONTROL IN AGED CARE

In accordance with Australian Aged Care accreditation standards¹, there are two key areas where Halyard products can help provide protection for residents, patients, staff and facilities:

- Patient Care and Contenance Management, and
- Infection Control.

Halyard's range of Clinical Towel, Fluid Absorption and Personal Protective products help to provide essential protection for Patients, Residents and Staff, as well as for surfaces and equipment in the Aged Care setting.

USE FOR:



PATIENTS



STAFF



SURFACES &
EQUIPMENT

HEALTH AND PERSONAL CARE:				
PATIENT CARE	HALYARD* Neutral Hand Wipes			page 4
	HALYARD* Protectapad			page 4
	HALYARD* Washcloth Extra			page 4
	HALYARD* Patient Mask			page 11
CONTINENCE MANAGEMENT	HALYARD* Underpads			page 5
	HALYARD* Bed Protector			page 5
INFECTION CONTROL:				
HAND HYGIENE	HALYARD* Neutral Hand Wipes			page 6
	HALYARD VERSA* Towel			page 7
	HALYARD TERI* Wipers			page 7
	HALYARD TERI* Hydro Absorbent Wipers			page 7
PERSONAL PROTECTIVE EQUIPMENT (PPE)	HALYARD* Masks			page 8
	HALYARD* Eyewear			page 12
	HALYARD* Gloves			page 14
	HALYARD* Gowns			page 20
	HALYARD* Footwear			page 24
SURFACE DISINFECTION	HALYARD ISOWIPE* Bactericidal Wipes			page 25
SURFACE PROTECTION	HALYARD* Underpads			page 26
	HALYARD* Bed Protector			page 26
	HALYARD* Bedsheet			page 26
	HALYARD* Benchroll			page 27
	HALYARD VERSA* Towel			page 27
	HALYARD TERI* Wipers			page 27
	HALYARD TERI* Hydro Absorbent Wipers			page 28
	HALYARD* Single Use Linen			page 29

PATIENT CARE

PATIENT HAND HYGIENE

HALYARD* NEUTRAL HAND WIPES



Alcohol-free wipes that contain neutral ingredients to clean visibly soiled hands when there is limited access to soap and water. They contain Vitamin E to restore moisture to skin for frequently washed hands.

- Ideal for limited mobility patients and their visitors
- Small soft pack is convenient, portable and easy to store
- Resealable pack helps reduce risk of contamination and drying out of wipes
- ARTG 273333

Code	Size	Description	Packaging
8835	20cm x 20cm	HALYARD* Neutral Hand Wipes	Eaches/Pack = 15 Packs/Case = 40 Eaches/Case = 600

PATIENT PERSONAL HYGIENE

HALYARD* PROTECTAPAD



Constructed from four highly absorbent layers of tissue and a layer of strong polyethylene backing, laminated together to help prevent strikethrough to provide absorbent waterproof surface protection. It is ideal for use as a bib to protect patients' clothing or for hygienic protection for patient surfaces.

- Embossed pattern increases liquid absorbency helping to prevent run off
- Laminated for extra strength, wet and dry²
- Softness ensures optimal patient comfort

Code	Size	Description	Packaging
2706	28.5cm x 21.5cm	HALYARD* Protectapad Small	Eaches/Pack = 100 Packs/Case = 8 Eaches/Case = 800
2705	28.5cm x 43cm	HALYARD* Protectapad Large	Eaches/Pack = 100 Packs/Case = 4 Eaches/Case = 400

HALYARD* WASHCLOTH EXTRA



An incredibly gentle, super absorbent washcloth that has been designed specifically for hygienic patient care. It's thick, soft and absorbs up to eight times its weight in water³. It is low-linting⁴ and is strong and durable enough to not tear easily² in wet or dry applications. HALYARD* Washcloth Extra is disposable, so it helps reduce the risks associated with cross-contamination, making it the hygienic choice for patient care and instrument protection.

- Ideal as face or body washer
- Clear resealable bag helps prevent dust and other particles from contaminating cloths

Code	Size	Description	Packaging
2711	35cm x 40cm	HALYARD* Washcloth Extra	Eaches/Pack = 20 Packs/Case = 30 Eaches/Case = 600

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

PATIENT CARE

PATIENT CONTINENCE MANAGEMENT



HALYARD* UNDERPADS

Constructed from a soft top layer and highly absorbent layers of tissue with a strong polyethylene backing that is quiet, drapable and helps prevent strikethrough of liquid. Available in two absorbency levels – 5 ply and 8 ply – so they provide hygienic protection for patients and surfaces across a large number of applications. Ideal for light to moderate fluid situations.

- Ideal for chair and mattress protection and as an underlay for medical procedures
- Soft, smooth poly base is quiet to assist with patient discretion and drapable to protect varied surfaces
- Blue folded edge provides a visual guide for edge of Underpad against bedsheets
- Hygienic protection that can assist in reducing the time taken for staff to clean and change linen

Code	Size	Description	Packaging
2765	39.5cm x 55cm	HALYARD* Underpads 5 Ply, Folded	Sheets/Polybag = 75 Polybags/Case = 4 Sheets/Case = 300
2768	39.5cm x 55cm	HALYARD* Underpads 8 Ply, Folded	Sheets/Polybag = 50 Polybags/Case = 4 Sheets/Case = 200



HALYARD* BED PROTECTOR

A unique design, constructed from four layers: a surface layer of soft spunbond, a second layer of special absorbent core which contains Super Absorbent Material, and a double layered base with strong poly-coated cover that helps prevent strikethrough. Ideal for moderate to heavy fluid situations.

- With a 4.9L absorbent capacity³, the highly absorbent core draws fluid away from the body and lock wetness away to help keep patients dry⁵ and reduce odour, providing ongoing absorbency
- Strong double layer poly base⁶ helps prevent strikethrough and maintain structural integrity when saturated, to keep surfaces dry and hygienic. Also helps reduce laundering costs and the need for staff to change bed linen
- Long ends⁶ provide extra material to secure pad in place

Code	Size	Description	Packaging
4261	81cm x 78cm	HALYARD* Bed Protector	Eaches/Case = 25

INFECTION CONTROL⁷

Infection prevention and control is an essential part of care and the responsibility of all staff providing care to older people. Infections can spread in any environment, however residents or clients in aged care may be more vulnerable to infection for a number of reasons including: being older, having been in hospital and having chronic diseases.

Communal living in residential care increases contact between people and also the risk of infection, especially when residents have difficulties with activities of daily living (e.g. toileting). In aged care facilities, residents and the staff members who care for them are the most likely sources of infectious agents to residents. In particular, care workers' hands are a likely source of transmission of infectious agents to residents.

Residents and care workers are the most common susceptible hosts in residential aged care facilities. Visitors and other people working in the facility may also be at risk of infection or transmission.

Care workers can help reduce the outbreak by rigorous adherence to hand hygiene, environmental cleaning protocols and appropriate use of Personal Protective Equipment (PPE).

HAND HYGIENE

Hand hygiene is the single most important factor in reducing the spread of infections. Hand hygiene can be performed either by washing hands with soap and water or using an alcohol-based hand rub.

In cases where access to soap and water is limited, healthcare workers and patients often turn to alcohol-based products to reduce the incidence of germs on their hands.

Alcohol is very effective in killing micro-organisms on the hands, however it's inactivated by organic matter and therefore not recommended as a cleaning agent⁸. Visibly soiled hands must be washed.



HALYARD* NEUTRAL HAND WIPES

Alcohol-free wipes that contain neutral ingredients to clean visibly soiled hands when there is limited access to soap and water. They contain Vitamin E to restore moisture to skin for frequently washed hands.

- Ideal for limited mobility patients and their visitors
- Small soft pack is convenient, portable and easy to store
- Resealable pack helps reduce risk of contamination and drying out of wipes
- ARTG 273333

Code	Size	Description	Packaging
8835	20cm x 20cm	HALYARD* Neutral Hand Wipes	Eaches/Pack = 15 Packs/Case = 40 Eaches/Case = 600

INFECTION CONTROL

HAND HYGIENE continued

HALYARD VERSA* TOWEL



VERSA* Towel is a high quality, white 2 ply crepe wadding towel that provides excellent absorbency and strength², so it can be used across a wide variety of medical and non-medical applications. As it's a soft and bulky crepe towel, it's gentle on skin without chafing from frequent hand drying. Ideal for use as soft absorbent hand towel or hygienic surface protection.

- Single-use towel helps prevent risk of cross contamination
- Wall or Bed Mounted Dispensers are available for easy and economical dispensing
- ARTG 332698

Code	Size	Description	Packaging
4210	24.5cm x 41.5cm	VERSA* Towel Small	Rolls/Case = 16
4220	49cm x 41.5cm	VERSA* Towel Large	Rolls/Case = 8
7041	27.8cm W x 14.8cm H x 14.5cm D	White Plastic Dispenser Suits Code 4210	1
7056	52.5cm W x 14.8cm H x 14.5cm D	White Plastic Dispenser Suits Code 4220	1

HALYARD TERI* WIPERS



TERI* Wipers are a soft 4 ply paper towel with a tough polyester scrim that provide excellent absorbency and incredible long-lasting strength even when wet². The thick 4 ply tissue provides a substantial barrier between your hands and contaminants and the soft embossed wiper is gentle on frequently washed hands.

- Single-use towel helps prevent risk of cross contamination
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4465	31.5 cm x 34cm	TERI* Wipers Small	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600
4467	31.5 cm x 34cm	TERI* Wipers Small Bulk Pack	Eaches/Pack = 450
4466	31.5 cm x 61cm	TERI* Wipers Large	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600

HALYARD TERI* HYDRO ABSORBENT WIPERS



TERI* Hydro Absorbent Wipers are made from highly absorbent material that can absorb up to six times its own weight in water³. The thick cloth-like feel is gentle on frequently washed hands and provides a substantial barrier between your hands and contaminants.

- Single-use – helps prevent risk of cross contamination
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4476	32.5cm x 33cm	TERI* Hydro Absorbent Wipers Small	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4478	32.5cm x 60cm	TERI* Hydro Absorbent Wipers Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

INFECTION CONTROL

PERSONAL PROTECTIVE EQUIPMENT (PPE)

Halyard offers products across all areas of Personal Protective Equipment (PPE) where there may be risk of exposure to blood or other bodily fluids.

MASKS⁷

Face masks are used to protect a care worker's nose and mouth from exposure to infectious agents.

They are used when there is a risk of:

- droplets or aerosols (e.g. from coughs or sneezes)
- splashes or sprays of blood or body fluids (e.g. when emptying wound or catheter bags).

Masks are worn during care of people who have an infection that is spread by the droplet or airborne route.

Masks may also be placed onto clients/residents who are coughing, especially if they are unable to cover their mouths. Before doing this, consider whether wearing a mask will cause distress (e.g. if the client or resident is unable to understand the purpose of wearing it).



GOLDEN RULES – MASKS:

- Check manufacturer's instructions before use.
- Don't touch front of the mask with your hands once the mask is in place.
- Use each mask for the care of one person only and change if a care activity is taking a long period of time.
- Don't leave mask dangling around your neck.
- Discard after use and perform hand hygiene after discarding.

HALYARD, YOUR PARTNER IN COMPLIANCE TO THE MASK STANDARD AS 4381:2015 FOR SINGLE-USE FACE MASKS

NEW STANDARD AS 4381:2015 (RED INDICATES CHANGE)				
AS 4381: 2015 SINGLE-USE FACE MASKS				
CHARACTERISTICS	LEVEL 1 1	LEVEL 2 2	LEVEL 3 3	TEST METHOD
	Level 1 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the minimum velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For general purpose medical procedures, where the wearer is not at risk of blood or bodily fluid splash or to protect staff and/or the patient from droplet exposure to microorganisms (e.g. patient with upper respiratory tract infection visits GP)	Level 2 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the middle velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For use in emergency departments, dentistry, changing dressings on small or healing wounds where minimal blood droplet exposure may possibly occur (e.g. endoscopy procedures)	Level 3 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the maximum velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For all surgical procedures, major trauma first aid or in any area where the health care worker is at risk of blood or bodily fluid splash (e.g. orthopaedic, cardiovascular procedures)	
Bacterial Filtration Efficiency (BFE) %	≥ 95%	≥ 98%	≥ 98%	ASTM F2101-14 or EN 14683:2014
Particulate Filtration Efficiency (PFE) % (0.1 µm)	Not Required	Not Required	Not Required	N/A
Differential Pressure (Delta P) mm H ₂ O/cm ²	< 4.0	< 5.0	< 5.0	EN 14683:2014
Resistance to penetration by synthetic blood (fluid resistance) min pressure in mm Hg for pass result	80mm Hg	120mm Hg	160mm Hg	ASTM F1862 / F1862M-13 or ISO 22609

INFECTION CONTROL

For the highest levels of protection, make sure you choose HALYARD* Face Masks and Eyewear.

NIOSH APPROVED N95 RESPIRATOR FACE MASKS

- Provides fluid resistance at 160mm Hg
- The highest level of protection available
- Recommended for Laser Procedures




46727

Code	Description	Packaging
46727	FLUIDSHIELD* Level 3 N95 Particulate Filter Respirator and Surgical Mask, Pouch-Style with Headbands, Regular Size, Orange	Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210
46827	FLUIDSHIELD* Level 3 N95 Particulate Filter Respirator and Surgical Mask, Pouch-Style with Headbands, Small Size, Orange	Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210

HALYARD FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 160mm Hg
- Level 3 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



28797



28800




28804




47107




47147



48207



48247



48297

Code	Description	Packaging
28797	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Ear Loops	Eaches/Box = 40 Boxes/Case = 10 Eaches/Case = 400
28800	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Wraparound Visor and SO SOFT* Lining and SO SOFT* Ear Loops	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
28804	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with Wraparound Visor and SO SOFT* Lining	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
47107	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask, Pleat-Style with SO SOFT* Ear Loops, Foam Band, Orange	Eaches/Box = 40 Boxes/Case = 10 Eaches/Case = 400
47147	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with SO SOFT* Ear Loops, Foam Band, Orange	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
48207	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask, Pleat-Style with Ties, Orange	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
48247	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties, Orange	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
48297	FLUIDSHIELD* Level 3 Care Bear Surgical Mask, Pleat-Style with Ties, Pink and Blue Teddy Bear Print	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

INFECTION CONTROL



HALYARD FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 120mm Hg
- Level 2 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



62114

Code	Description	Packaging
62114	FLUIDSHIELD* Level 2 Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties, Blue with Orange Diamond Print	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100



62115

62115	FLUIDSHIELD* Level 2 Fog-Free Procedure Mask with SO SOFT Ear Loops, Blue	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
-------	---	--



HALYARD FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 80mm Hg
- Level 1 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



25867

Code	Description	Packaging
25867	FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining and Earloops, Yellow	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500



25868

25868	FLUIDSHIELD* Level 1 Procedure Mask, Pleat-Style with SO SOFT* Ear Loops, Lavender	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	--	---



28802

28802	FLUIDSHIELD* Level 1 Surgical Mask with SO SOFT* Lining	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	---	---



41802

41802	FLUIDSHIELD* Level 1 Procedure Mask with Fog-Free Foam Strip and SO SOFT* Lining and Ear Loops, White	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	---	---

INFECTION CONTROL

STANDARD FACE MASKS



HALYARD* PROCEDURE MASKS

- Protective three-layer mask construction, for when exposure to blood and/or bodily fluids is not a risk
- Appropriate for general patient care and isolation
- Earloop attachment

Code	Description	Packaging
47080	Standard Procedure Mask, Earloops, Blue	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500



MASK DISPENSERS

HALYARD* MASK DISPENSER WALL BRACKET

Code	Description	Packaging
4966	Wall Bracket, 170mm (L) x 100mm (W) x 160mm (H)	Eaches = 1
4967	Wall Bracket, 170mm (L) x 125mm (W) x 160mm (H)	Eaches = 1
4968	Wall Bracket, 170mm (L) x 190mm (W) x 160mm (H)	Eaches = 1

INFECTION CONTROL


PROTECTIVE EYEWEAR⁷

Protective eyewear is used to protect a care worker's eyes from exposure to infectious agents.

It is used when there is a risk of:

- droplets or aerosols (e.g. from oral suctioning)
- splashes or sprays of blood or body fluids (e.g. when emptying catheter bags).

Eyewear should be worn during care of people who have an infection that is spread by the droplet or airborne route.



GOLDEN RULES – PROTECTIVE EYEWEAR:

- Remember that the outside of the eyewear is contaminated.
- Remove using the headband or ear pieces.
- Clean eye shield after each use with detergent and water and allow it to dry.
- If eyewear is single use, dispose after the care activity.

EYEWEAR



HALYARD GUARDALL* FACE SHIELD

- Designed for procedures and environments where risk of exposure to blood and/or bodily fluids through splashing is possible
- Equipped with a foam band and headband
- Can be worn with or without glasses or goggles
- Resistant to fogging

Code	Description	Packaging
41204	GUARDALL SHIELD* Full-Length Visor Only, (No Mask), Clear	Eaches/Case = 40
41205	GUARDALL SHIELD* Three-Quarter Length Visor Only (No Mask), Clear	Eaches/Case = 40



HALYARD SAFEVIEW* GLASSES

- For use in situations which may involve splashing or spraying
- Designed to be lightweight for added comfort in use
- Shipped with colorful frames to appeal to a variety of preferences
- Single-use lenses may be easily replaced into reusable frames
- Designed to be easily worn over most prescription glasses
- Constructed from natural rubber latex-free materials



Code	Description	Packaging
SV50A	SAFEVIEW* Assembled Glasses	5 Pkg. of 10 Each Case = 50 Assembled Glasses
SV100F	SAFEVIEW* Frames / Assorted Colours	10 Pkg. of 10 Each Case = 100 Frames
SV250L	SAFEVIEW* Replacement Lenses	10 Pkg. of 25 Each Case = 250 Lenses

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

INFECTION CONTROL

TEST METHODS FOR FACE MASKS

BACTERIAL FILTRATION EFFICIENCY (BFE) AT 3 μ m: ASTM F2101

Determines the bacterial filtration efficiency (BFE) of various filtration materials. This test employs a ratio of the bacterial challenge counts to sample effluent counts to determine the percent bacterial filtration efficiency (%BFE). The test procedure allows a reproducible bacterial challenge to be delivered to test materials, and provides a standard procedure for comparison of filtration materials.

The higher the BFE, the better the filtration efficiency, the better the protection.

DIFFERENTIAL PRESSURE OR BREATHABILITY (DELTA P): 4.4.1.2 of MIL-M-36954C

Determines the air exchange differential of porous materials measured in mm H₂O/cm². This technique employs a water manometer differential upstream and downstream of the test material at a constant flow rate.

The lower the delta P, the more breathable the material is. However, this can also mean lower levels of protection.

FLUID RESISTANCE: ASTM F 1862

Standard Test Method for Resistance of Medical Face Masks to Penetration by Synthetic Blood (Horizontal Projection of Fixed Volume at Known Velocity).

Developed to simulate a challenge to the fluid splash resistance of a face mask under conditions similar to actual use. During this test, face masks are insulted with synthetic blood, and graded as a pass or fail. The masks are evaluated at three pressures: 80, 120, and 160mmHg.

The higher the pressure at which a mask passes, the greater the fluid resistance.

NATIONAL INSTITUTE FOR OCCUPATIONAL SAFETY AND HEALTH (NIOSH) Certification for Particulate Respirators: 42 CFR 84 Certified in USA

NIOSH Certification ensures the respirator has undergone the stringent testing and quality assurance requirements to meet an expected level of protection.

If a respirator has been NIOSH-certified, it will have the following markings on the respirator's packaging, user instruction insert, and/or on the respirator itself:

- The NIOSH name or logo;
- The NIOSH filter classification – in this case "N95";
- The NIOSH test and certification – or TC – approval number;
- The Approval Holder's business name (usually the manufacturer), registered trademark, or an easily understood abbreviation;
- The model or part number; and
- The manufacturing lot number.

The respirators certified under the NIOSH 42 CFR 84 regulation are tested under much more demanding conditions, providing a higher level of worker protection.

INFECTION CONTROL

GLOVES⁷

Gloves are used to reduce the spread of infectious agents that may be carried on hands.

Gloves are worn for procedures where there is a risk of exposure to blood or body substances.

They are also worn during care of people who have an infection spread by contact (especially if it is caused by a multi-resistant organism).



Wear gloves when:

- **Changing a colostomy bag or urinary drainage bag**
- **Dressing wounds or touching broken skin**
- **Assisting with toileting**
- **Giving mouth or eye care**
- **Oral suctioning**
- **Touching equipment or surfaces that may come into contact with blood or body substances**
- **Blood glucose monitoring**
- **You have broken skin**
- **Preparing food.**

GOLDEN RULES – GLOVES:

- **Gloves are not used instead of hand hygiene.**
- **Perform hand hygiene before and after using gloves.**
- **Remove gloves when a care activity is finished. Change gloves before starting a different care activity.**
- **Dispose of used gloves immediately.**

INFECTION CONTROL

THE WORRY-FREE CHOICE FOR EVERYONE, EVERYWHERE

Your hands are essential to the care you provide to your patients. With 5–15% of healthcare workers susceptible to latex allergy, caring for you and your hands by reducing the exposure to latex reduces the risks and harmful effects of an allergic reaction.⁹

Halyard has created a patented, natural rubber latex-free glove to provide you with the comfort, fit and sensitivity of latex, without the added risk of allergy exposure.

It's the one glove for all your non-surgical needs, facility-wide.¹⁰



INNOVATION IN PACKAGING TO REDUCE CROSS CONTAMINATION AND WASTE

As part of our commitment to reducing the risk of HAIs, Halyard introduced the SmartPULL* glove dispensing system to our STERLING* and LAVENDER* Nitrile Exam Glove boxes.

Our advanced packing technology has also reduced the waste in storage space with the larger, 200–250 count per box providing you with more for less.



The smaller tab reduces the risk of contamination and waste by limiting the contact and potential fall-out from adjacent gloves.



The larger tab enables easier access as the box empties.

INFECTION CONTROL

NITRILE EXAM GLOVES



HALYARD STERLING* NITRILE EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- Silver
- CE Marked
- Enhanced Wet and Dry Grip
- Tested Against 12 Chemotherapy Agents



Code	Description	Packaging
13938	STERLING* Nitrile Exam Gloves, X-Small	Eaches/Box = 200 Boxes/Case = 10
13940	STERLING* Nitrile Exam Gloves, Small	Eaches/Box = 200 Boxes/Case = 10
13941	STERLING* Nitrile Exam Gloves, Medium	Eaches/Box = 200 Boxes/Case = 10
13942	STERLING* Nitrile Exam Gloves, Large	Eaches/Box = 200 Boxes/Case = 10
13943	STERLING* Nitrile Exam Gloves, X-Large	Eaches/Box = 170 Boxes/Case = 10

HALYARD STERLING* NITRILE-XTRA EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 12" (305mm)
- Silver
- Enhanced Wet and Dry Grip
- Tested Against 12 Chemotherapy Agents



Code	Description	Packaging
53137	STERLING* NITRILE-XTRA Exam Gloves, X-Small	Eaches/Box = 100 Boxes/Case = 10
53138	STERLING* NITRILE-XTRA Exam Gloves, Small	Eaches/Box = 100 Boxes/Case = 10
53139	STERLING* NITRILE-XTRA Exam Gloves, Medium	Eaches/Box = 100 Boxes/Case = 10
53140	STERLING* NITRILE-XTRA Exam Gloves, Large	Eaches/Box = 100 Boxes/Case = 10
53141	STERLING* NITRILE-XTRA Exam Gloves, X-Large	Eaches/Box = 100 Boxes/Case = 10

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

INFECTION CONTROL



HALYARD PURPLE NITRILE* EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- Purple
- CE Marked
- Enhanced Wet and Dry Grip
- Tested Against 29 Chemotherapy and other Hazardous Agents



Code	Description	Packaging
52000	PURPLE NITRILE* Exam Gloves, X-Small	Eaches/Box = 100 Boxes/Case = 10
52001	PURPLE NITRILE* Exam Gloves, Small	Eaches/Box = 100 Boxes/Case = 10
52002	PURPLE NITRILE* Exam Gloves, Medium	Eaches/Box = 100 Boxes/Case = 10
52003	PURPLE NITRILE* Exam Gloves, Large	Eaches/Box = 100 Boxes/Case = 10
52004	PURPLE NITRILE* Exam Gloves, X-Large	Eaches/Box = 90 Boxes/Case = 10



HALYARD PURPLE NITRILE-XTRA* HIGH RISK EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 12" (305mm)
- Longer Cuff for Extra Protection
- CE Marked
- Purple
- Enhanced Wet and Dry Grip
- Tested Against 29 Chemotherapy and other Hazardous Agents



Code	Description	Packaging
50600	PURPLE NITRILE-XTRA* Exam Gloves, X-Small	Eaches/Box = 50 Boxes/Case = 10
50601	PURPLE NITRILE-XTRA* Exam Gloves, Small	Eaches/Box = 50 Boxes/Case = 10
50602	PURPLE NITRILE-XTRA* Exam Gloves, Medium	Eaches/Box = 50 Boxes/Case = 10
50603	PURPLE NITRILE-XTRA* Exam Gloves, Large	Eaches/Box = 50 Boxes/Case = 10
50604	PURPLE NITRILE-XTRA* Exam Gloves, X-Large	Eaches/Box = 50 Boxes/Case = 10

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

INFECTION CONTROL



HALYARD PURPLE NITRILE* MAX EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Purple
- Longer length: 16" (406mm)
- CE Marked
- Enhanced Wet and Dry Grip
- Heavier weight: 20 grams
- Increased fingertip thickness
- Tested Against 52 Chemotherapy and other Hazardous Agents



Code	Description	Packaging
44992	PURPLE NITRILE* MAX Exam Gloves, Small	Eaches/Box = 50 Boxes/Case = 8
44993	PURPLE NITRILE* MAX Exam Gloves, Medium	Eaches/Box = 50 Boxes/Case = 8
44994	PURPLE NITRILE* MAX Exam Gloves, Large	Eaches/Box = 50 Boxes/Case = 8
44995	PURPLE NITRILE* MAX Exam Gloves, X-Large	Eaches/Box = 50 Boxes/Case = 8



HALYARD LAVENDER* NITRILE EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- Lavender
- Enhanced Wet and Dry Grip

Code	Description	Packaging
52816	LAVENDER* Nitrile Exam Gloves, X-Small	Eaches/Box = 250 Boxes/Case = 10
52817	LAVENDER* Nitrile Exam Gloves, Small	Eaches/Box = 250 Boxes/Case = 10
52818	LAVENDER* Nitrile Exam Gloves, Medium	Eaches/Box = 250 Boxes/Case = 10
52819	LAVENDER* Nitrile Exam Gloves, Large	Eaches/Box = 250 Boxes/Case = 10
52820	LAVENDER* Nitrile Exam Gloves, X-Large	Eaches/Box = 230 Boxes/Case = 10

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

INFECTION CONTROL



50201



50202



50203

EVERYDAY GLOVE DISPENSERS

Code	Description	Packaging
50201	Exam Glove Dispenser, Clear Plastic, Single, 96mm (D) x 144mm (W) x 253mm (H)	Eaches = 1
50202	Exam Glove Dispenser, Clear Plastic, Double, 96mm (D) x 274mm (W) x 256mm (H)	Eaches = 1
50203	Exam Glove Dispenser, Clear Plastic, Triple, 96mm (D) x 406mm (W) x 256mm (H)	Eaches = 1



47467



47466

HALYARD* PURPLE NITRILE* MAX GLOVE DISPENSERS

Code	Description	Packaging
47467	Exam Glove Dispenser Extra Long (Suits 16" gloves),, Clear Plastic, Single, 75mm (D) x 139mm (W) x 409mm (H)	Eaches = 1
47466	Exam Glove Dispenser Extra Long (Suits 16" gloves),, Clear Plastic, Triple, 69 mm (D) x 399mm (W) x 409 mm (H)	Eaches = 1


INFECTION CONTROL

GOWNS⁷

Impermeable (waterproof) gowns or aprons are used to stop contamination of care workers' clothes and skin.

They are used when there is a risk of splashes or sprays of blood or body fluids (e.g. if there is vomiting or diarrhoea).

Gowns/aprons are worn during care of people who have an infection that is spread by the contact, droplet or airborne route.



GOLDEN RULES – GOWNS:

- Hand hygiene must be performed before and after using gowns or aprons.
- Gowns must fully cover the torso from neck to knees, arms to end of wrists and adequately wrap around the back. Tie all fastenings on the gown and fasten at the back.
- Remove and dispose of the gown as soon as care is completed.

SEQUENCE FOR **DONNING** PPE:

1. PROTECT YOUR BODY

Ensure the gown:

- Covers torso from neck to knees
- Covers arms to end of wrists
- Wraps around the sides to fasten at neck and back



2. PROTECT YOUR FACE

- Create a curve in the flexible band for the nose bridge and expand the mask
- Secure the ties or elastics at the middle of the head and neck, or fasten the ear loops
- Fit it snugly to your face and below your chin – fit check the N95 and P2 mask
- Place the face shield or eyewear over the face and eyes and adjust to fit



3. PROTECT YOUR HANDS

- Extend the cuffs to cover the wrist of the gown
- Limit the surfaces touched
- Change gloves when torn or heavily contaminated
- Perform hand hygiene as required



FINALLY STOP and CHECK all PPE is in place



INFECTION CONTROL

SEQUENCE FOR **ADVANCED DOFFING PPE:**

Remember! The outside and front of your PPE, including the sleeves of the gown, are contaminated. If hands become visibly contaminated during PPE removal, wash hands before continuing.

REMOVE **GOWN THEN GLOVES**

Outside of gown and gloves are CONTAMINATED!

1a.

- With gloved hands, grasp the front of the gown at the shoulders and pull the gown forward to snap the tabs at the back
- Peel off the gown, rolling it inside out

1b.

- Continue over the gloves so that they too are rolled inside out
- Discard gown and gloves into a waste container
- Perform hand hygiene using ABHR or soap and water



2. REMOVE **PROTECTIVE EYEWEAR**

- Grasp the ear or head piece and lift away from face
- Place reusable eyewear in the designated receptacle for cleaning or
- Discard disposable eyewear into waste container for disposal



3. REMOVE **MASK OR PARTICULATE RESPIRATOR (N95 OR P2)**

- Remove by touching ties or elastics only, starting with the bottom then top ties or bands
- Discard in designated waste container



4. **PERFORM HAND HYGIENE**

- Wash hands with soap and water or use an alcohol-based hand rub immediately after removing PPE



INFECTION CONTROL

HALYARD* ISOLATION GOWNS

- For use when expected risk of exposure to fluid is between low and moderate
- Medium-weight tri-layer fabric
- Tape-tab neck closure
- Tie waist, elastic cuffs



Code	Description	Packaging
69979	Isolation Gown, Yellow, Large	Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100
69988	Isolation Gown, Yellow, X-Large	Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100



HALYARD* IMPERVIOUS COMFORT GOWNS

- Potential for moderate to heavy fluid contact in front
- Tri-layer spunbond with polyethylene laminate coating
- Breathable spunbond back
- Knit cuffs
- Velcro neck tabs
- Tie waist



Code	Description	Packaging
69455	Impervious Open Back Gown, Knit Cuff, Blue, Large	Eaches/Bundle = 30 Bundles/Case = 2 Eaches/Case = 60
5511	Impervious Full Back Comfort Gown, Yellow, Regular	Eaches/Case = 50
5512	Impervious Full Back Comfort Gown, Yellow, X-Large	Eaches/Case = 50

INFECTION CONTROL



HALYARD* IMPERVIOUS THUMBS-UP* FILM GOWNS

- Potential for soaking and heavy fluid contact in front
- Impervious polyethylene film
- Open back, thumbhooks

Code	Description	Packaging
7000	Impervious Thumbs-Up* Film Gown with Thumbhooks, Blue, Regular	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7001	Impervious Thumbs-Up* Film Gown with Thumbhooks, Blue, X-Large	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7002	Impervious Thumbs-Up* Film Gown with Thumbhooks, Yellow, Regular	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7003	Impervious Thumbs-Up* Film Gown with Thumbhooks, Yellow, X-Large	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75



HALYARD* CHEMOTHERAPY GOWN

- Potential for soaking and heavy fluid contact in front
- Tested with the 7 chemotherapy drugs required by ASTM test method F739
- Tested with 3 additional standard widely used chemotherapy drugs: meets Liquid Barrier Performance Standards ASTM 1670 and 1671
- Highest level of protection against flammability
- Polypropylene with polyethylene coating



Code	Description	Packaging
69606	Comfort Gown, Knit Cuffs, Perforated Open Back, Blue, Universal	Eaches/Case = 100
37284	Comfort Gown, Knit Cuffs, Perforated Open Back, Blue, XX-Large	Eaches/Case = 100



HALYARD* CHEMO360* PROCEDURE GOWN

- Closed in the back
- Secure tape closure

Code	Description	Packaging
47346	CHEMO360* Gown, Closed Back, Poly-Coated, Tape Closure, Blue, Universal	Eaches/Case = 100
47347	CHEMO360* Gown, Closed Back, Poly-Coated, Tape Closure, Blue, XX-Large	Eaches/Case = 100

INFECTION CONTROL

FOOTWEAR



HALYARD* OVERSHOES

Code	Description	Packaging
69252	X-TRA TRACTION* Shoe Cover, Universal	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300
69254	X-TRA TRACTION* Shoe Cover, X-Large E	Eaches/Dispenser = 80 Dispensers/Case = 3 Eaches/Case = 240

ANKLE TOPS



HALYARD ANKLE GUARD* SHOE COVER

- Moderate fluid contamination
- Over the ankle coverage
- Protective tri-layer SMS fabric
- Adhesive traction strips
- Natural rubber latex-free

Code	Description	Packaging
69353	ANKLE GUARD* Shoe Cover, Tri-Layer Fabric, Adhesive Traction Strips, X-Large	Eaches/Dispenser = 90 Dispensers/Case = 3 Eaches/Case = 270

KNEE HIGHS



HALYARD HI-GUARD* ULTRA FULL-COVERAGE BOOT

- Potential for heavy fluid contamination
- Protection for staff when assisting with showering
- Knee high boot
- Plastic film coating over tri-layer fabric
- Two 1", 2.5cm adhesive foam traction strips
- Natural rubber latex-free

Code	Description	Packaging
69572	HI-GUARD* Full Coverage Boot, Medium	Eaches/Dispenser = 30 Dispensers/Case = 4 Eaches/Case = 120
69672	HI-GUARD* Full Coverage Boot, Large	Eaches/Dispenser = 30 Dispensers/Case = 4 Eaches/Case = 120

INFECTION CONTROL

PROTECTION FOR SURFACES AND EQUIPMENT

CLEANING AND SURFACE DISINFECTION⁷

The level of cleaning required to stop the spread of infection depends on the objects involved and the risk of contamination.

Extra measures (e.g. wearing gloves and gowns, cleaning and disinfecting surfaces likely to be touched by the resident/client) may be needed if the person has an uncovered wound or has trouble with hand hygiene or cough etiquette.

GOLDEN RULES – CLEANING

- Most surfaces can be adequately cleaned with warm water and detergent per manufacturer's instructions.
- Allow the cleaned surface to dry completely.
- Detergent solution followed by disinfectant may be appropriate when an infection is known or suspected.

SURFACE DISINFECTION

HALYARD ISOWIPE* BACTERICIDAL WIPES



Containing 70% Isopropyl Alcohol, ISOWIPE* Bactericidal Wipes are a TGA-listed hospital grade disinfectant that offers disinfection of clean surfaces, in a convenient wet wipe that doesn't require dilution. The foil seal helps to maintain the integrity of the formula during storage and transportation and the specially designed dispensing slot controls dispensing and minimises the evaporative effect of the alcohol.

- Contains 70% Isopropyl Alcohol to provide hospital grade disinfection of clean surfaces
- Kills Bacteria (*S. aureus*, *E. coli*, *Ps. Aeruginosa*, *Pr. Vulgaris*)
- Kills Fungi (Yeast and Mould) in 10 minutes
- HACCP certified
- Entered in the ARTG as hospital grade disinfectant ARTG 342811

Always read the label and follow the directions for use.



Code	Size	Description	Packaging
6835	42cm x 14cm	ISOWIPE* Bactericidal Wipes Canister	Eaches/Canister = 75 Canisters/Case = 12 Eaches/Case = 900
6836	42cm x 14cm	ISOWIPE* Refill Pack	Eaches/Refill = 75 Refills/Case = 12 Eaches/Case = 900
6837	21cm x 14cm	ISOWIPE* Minis Canister	Eaches/Canister = 150 Canisters/Case = 12 Eaches/Case = 1800
6838	21cm x 14cm	ISOWIPE* Minis Refill Pack	Eaches/Refill = 150 Refills/Case = 12 Eaches/Case = 1800
4989	17.8cm H x 9.4cm D	Stainless Steel Wall Bracket	1

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

INFECTION CONTROL

SURFACE PROTECTION



HALYARD* UNDERPADS

Constructed from a soft top layer and highly absorbent layers of tissue with a strong polyethylene backing that is quiet, drapable and helps prevent strikethrough of liquid. Available in two absorbency levels – 5 ply and 8 ply – so they provide hygienic protection for patients and surfaces across a large number of applications. Ideal for light to moderate fluid situations.

- Ideal for chair and mattress protection and as an underlay for medical procedures
- Soft, smooth poly base is quiet to assist with patient discretion and drapable to protect varied surfaces
- Blue folded edge provides a visual guide for edge of Underpad against bedsheet
- Hygienic protection that can assist in reducing the time taken for staff to clean and change linen

Code	Size	Description	Packaging
2765	39.5cm x 55cm	HALYARD* Underpads 5 Ply, Folded	Eaches/Case = 300
2768	39.5cm x 55cm	HALYARD* Underpads 8 Ply, Folded	Eaches/Case = 200



HALYARD* BED PROTECTOR

A unique design, constructed from four layers: a surface layer of soft spunbond, a second layer of special absorbent core which contains Super Absorbent Material, and a double layered base with strong poly-coated cover that helps prevent strikethrough. Ideal for moderate to heavy fluid situations.

- With a 4.9L absorbent capacity³, the highly absorbent core draws fluid away from the body and lock wetness away to help keep patients dry⁵ and reduce odour, providing ongoing absorbency
- Strong double layer poly base⁶ helps prevent strikethrough and maintain structural integrity when saturated, to keep surfaces dry and hygienic. Also helps reduce laundering costs and the need for staff to change bed linen
- Long ends⁶ provide extra material to secure pad in place

Code	Size	Description	Packaging
4261	81cm x 78cm	HALYARD* Bed Protector	Eaches/Case = 25



HALYARD* BEDSHEET

Bedsheet is a thick durable tissue with a linen-like feel providing optimal patient comfort and hygienic surface protection for examination beds. The continuous roll can be cut to desired length and is disposable, which saves you time and effort in laundering and keeps your beds looking as good as new.

Code	Size	Description	Packaging
4260	54cm x 80m	HALYARD* Bedsheet Roll	Rolls/Case = 6
4924		Dispenser White Steel Suits Code 4260	1

INFECTION CONTROL



HALYARD* BENCHROLL

Constructed from four highly absorbent layers of tissue and one layer of strong polyethylene backing, laminated together for extra strength² wet or dry to help prevent strikethrough. It is the ideal choice for hygienic protection by protecting bench and surfaces from staining and other spills.

- Embossed pattern increases liquid absorbency so spills are quickly soaked up, helping to prevent run off
- Disposable – helps prevent risk of cross contamination
- Economical – cut to desired length to reduce unnecessary waste



Code	Size	Description	Packaging
2911	41.5cm x 91m	HALYARD* Benchroll 4 Ply Roll	Rolls/Case = 2



HALYARD VERSA* TOWEL

VERSA* Towel is a high quality, white 2 ply crepe wadding towel that provides excellent absorbency and strength², so it can be used across a wide variety of medical and non-medical applications. As it's a soft and bulky crepe towel, it's gentle on skin without chafing from frequent hand drying. Ideal for use as soft absorbent hand towel or hygienic surface protection.

- Single-use towel helps prevent risk of cross contamination
- Wall or Bed Mounted Dispensers are available for easy and economical dispensing
- ARTG 332698



Code	Size	Description	Packaging
4210	24.5cm x 41.5cm	VERSA* Towel Small	Rolls/Case = 16
4220	49cm x 41.5cm	VERSA* Towel Large	Rolls/Case = 8
7041	27.8cm W x 14.8cm H x 14.5cm D	White Plastic Dispenser Suits Code 4210	1
7056	52.5cm W x 14.8cm H x 14.5cm D	White Plastic Dispenser Suits Code 4220	1



HALYARD TERI* WIPERS

TERI* Wipers are a soft 4 ply paper towel with a tough polyester scrim that provide excellent absorbency and incredible long-lasting strength even when wet². The thick 4 ply tissue provides a substantial barrier between your hands and contaminants and the soft embossed wiper is gentle on frequently washed hands.

- Single-use towel helps prevent risk of cross contamination
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698



Code	Size	Description	Packaging
4465	31.5 cm x 34cm	TERI* Wipers Small	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600
4467	31.5 cm x 34cm	TERI* Wipers Small Bulk Pack	Eaches/Pack = 450
4466	31.5 cm x 61cm	TERI* Wipers Large	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600

INFECTION CONTROL



HALYARD TERI* HYDRO ABSORBENT WIPERS

TERI* Hydro Absorbent Wipers are made from highly absorbent material that can absorb up to six times its own weight in water³. The thick cloth-like feel is gentle on frequently washed hands and provides a substantial barrier between your hands and contaminants.

- Single-use – helps prevent risk of cross contamination
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4476	32.5cm x 33cm	TERI* Hydro Absorbent Wipers Small	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4478	32.5cm x 60cm	TERI* Hydro Absorbent Wipers Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450

PATIENT LINEN



HALYARD* SINGLE USE LINEN

Code	Description	Size	Packaging
67783	Stretcher Sheet, Blue	101.5cm x 213cm (40" x 84"), Bulk	80 eaches/case
67843	Draw Sheet, Poly-Coated*, Blue <small>*Coated with Polyethylene Film</small>	101.5cm x 183cm (40" x 72"), Bulk	80 eaches/case
67803	Pillowcase, Blue	50.5cm x 73.5cm (20in x 29in), Bulk	200 eaches/case

**For more information, please contact your local
O&M Halyard Account Manager or Customer Service on:**

 1800 664 227

 CustomerService_ANZ@hyh.com

 <https://www.halyardhealth.com.au/products/infection-prevention/aged-care.aspx>

References

1. Quality of Care Principles 2014, Aged Care Act, Australian Government.
2. Mechanical break resistance test data on file.
3. Absorption capacity test data on file.
4. Gelbo lint test data on file.
5. Re-wet test data on file. Note: These tests are directional based on small sample sizes.
6. Halyard Product Specifications data on file.
7. Australian Guidelines for the Prevention and Control of Infection in Healthcare, Canberra: National Health and Medical Research Council (2019).
8. Russell, A.D., Hugo, W.B., Ayliffe, G.A.J., Principles and Practice of Disinfection, Preservation & Sterilization, Fourth Edition, Blackwell Scientific Publication, London, 2004.
9. 2019, The Australasian Society of Clinical Immunology and Allergy (ASCIA) https://www.allergy.org.au/images/pcc/ASCIA_PCC_Latex_allergy_2019.pdf
10. Halyard offers a line of gloves specifically for use with chemotherapeutic agents.

O&M Halyard makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2022 O&M Halyard, Inc. All Rights Reserved.